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1008274001-034J  
Rev. 03/24AIA-PACK HBV Antibody  
CONTROL SET

## INTENDED USE

The AIA-PACK HBV Antibody CONTROL SET is intended for IN VITRO DIAGNOSTIC USE ONLY for performing quality control procedures with the ST AIA-PACK HBV antibodies assay.

## CONTENTS

Cat. No. 0020950  
2 x 4 mL AIA-PACK HBV Antibody CONTROL LEVEL 1  
2 x 4 mL AIA-PACK HBV Antibody CONTROL LEVEL 2

## WARNING AND PRECAUTIONS

1. The AIA-PACK HBV Antibody CONTROL SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. This material contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such materials, always flush with large volumes of water to prevent azide build up.
4. Although material derived from human blood is not used for the control set, it is recommended that this product be handled with the same precautions as used for patient samples.
5. Do not use beyond the expiration date.
6. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

## PREPARATION OF REAGENTS

The AIA-PACK HBV Antibody CONTROL SET is provided ready for use. Bring the controls to 18-25℃ prior to use.

## STORAGE AND STABILITY

Always store the AIA-PACK HBV Antibody CONTROL SET in an upright position at 2-8℃ when not in use. When stored unopened and refrigerated at 2-8℃, the control set is stable until the expiration date on the label. Once opened, the control materials should be used within 7 days provided that the vials are kept tightly sealed and refrigerated at 2-8℃.

## PROCEDURE


Set up the TOSOH AIA System Analyzers for analysis. Refer to the TOSOH AIA System Operator's Manual.

## EXPECTED VALUES AND LIMITATIONS

The ranges in the Data Sheet on the reverse side were obtained by using ST AIA-PACK available at the time of testing. The ranges are served as guidelines only for this specified lot of AIA-PACK HBV Antibody CONTROL SET, since the obtained target ranges were estimated by limited data. Since the assay values are dependent upon assay procedures as well as several other factors, each laboratory should establish its own range for the assay procedure being monitored.

TOSOH CORPORATION  
2-2-1, Yaesu, Chuo-ku,  
Tokyo 104-0028 Japan  
Phone: +81 3 6636 3734  
Fax : +81 3 6636 3627

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