

東ソーコントロールセット（心疾患）参考値一覧

2019.11.15 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)		CL-cTnI (pg/mL)		CL-CKMB (ng/mL)		CL-MYO (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
J249651	29-44	400-590	0.35-0.52	4.0-6.0	7.3-11	79-120	44-66	230-360	25-38	330-490	260-390	3200-4800	6.8-10	71-110	37-55	200-300
J649652	30-45	380-570	0.36-0.54	3.9-5.9	7.7-12	82-120	41-61	230-350	26-39	330-490	240-350	2900-4400	7.1-11	74-110	34-51	210-310
JX49653	29-43	390-580	0.37-0.56	3.9-5.9	7.9-12	87-130	41-62	220-330	27-41	340-510	260-390	3100-4600	7.3-11	77-110	36-55	200-300

Lot.	使用期限
J249651	2020.02
J649652	2020.06
JX49653	2020.10